


## S-CAPE® viewing solution

- S-CAPE® dicom viewer
- S-CAPE® collector
- S-CAPE® lizenzmanger



### Software-Packet for all complex solutions for

**S-CAPE®/PLANILUX® medical multiconsole**  
**S-CAPE®/PLANILUX® viewing station**

- Image viewing system for the OR
- Equipment for image viewing

Software-package      **S-CAPE® viewing solution**      enlarged the ambit of  
                                  **S-CAPE® dicom viewer**      via connection of  
                                  **S-CAPE® collector**      and  
                                  **S-CAPE® lizenzmanger**      to a comfortable multi user system.

The main part is      **S-CAPE® dicom viewer**

Software, which provides that all DICOM patients CD's could be evaluate. No matter of any difference from which CD manufacture the patient information comes from. The implementation of the S-CAPE® dicom viewer will provide a time saving of app. 40%.

Often users like to have a temporary limited saving of image data.

#### **S-CAPE® collector**

Enable this for non-medical workplaces (like the front desk).

With the S-CAPE® collector less time for requesting the patient's data is needed; it will reduce the time factor by 2 to 5%.

For cost reductions:

#### **S-CAPE® license-manager**

Is used for controlling the handling of the various Licenses (concurrent client) in use.

The installation of S-CAPE® viewer Software should be installed at all essential workstations.

Last Revision October 2008

